

# DNA-Miniprep. - Rapid boiling

by A. Untergasser (contact address and download at [www.untergasser.de/lab](http://www.untergasser.de/lab))  
Version: 1.0 - [Print Version \(.PDF\)](#)

**ATTENTION:** This is a low priced protocol. Use it preferably!

1. Pick colony and grow in 3 ml LB + antibiotic over night at 37°C
2. Prepare the STET++ and switch on the heat block to 95 °C
3. Spin down 1.5 ml for **2 min at 8000 rpm** (6000 G) store the rest at 4 °C
4. Resuspend pellet in **300 µl STET++** buffer
5. Incubate at room temperature for **10 min**
6. Incubate at **95 °C** (in heat block) for **2 min (not longer!)**
7. Spin down for **15 min** at 13000 rpm (18000 G)
8. **Discard pellet(!)** with a toothpick a toothpick
9. Add **500 µl Isoprop/NH4Ac** and mix
10. **Incubate for 20 min** at room temperature
11. Spin down for **15 min** at 13000 rpm (18000 G)
12. Discard supernatant and wash with **500 µl 70% Ethanol**
13. Spin down for **10 min** at 13000 rpm (18000 G)
14. Discard supernatant
15. Spin down for **2 min** at 13000 rpm (18000 G)
16. Pipett remaining supernatant of
17. Let pellet dry for 1 min
18. Resuspend pellet in **30 µl demi water**
19. **Digest 5 µl** of DNA solution in a total volume of 30 µl (**MasterMix !**)
20. Load **10 µl** of the digest on agarose gel

## **Buffers:**

### **STET:**

8 % Sucrose (w/v)  
0.1 % Triton X-100  
50 mM EDTA pH 8.0  
50 mM TrisCl pH 8.0

### **Isoprop/NH4Ac:**

75 ml Isopropanol 100%  
19.3 g Ammoniumacetat (NH4Ac)

### **LB:**

10 g Tryptone  
5 g Yeast Extract  
5 g NaCl

add water to **100 ml**  
(takes some time to dissolve)

add water to **1 litre**

## **Stock solutions:**

Lysozyme 50 mg / ml  
RNase A 10 mg / ml

## **Prepare STET++:**

	1x	10x	15x	20x	30x	60x	90x
STET	300 µl	3 ml	4.5 ml	6 ml	9 ml	18 ml	27 ml
Lysozym	5 µl	50 µl	75 µl	100 µl	150 µl	300 µl	450 µl
RNase A	1 µl	10 µl	15 µl	20 µl	30 µl	60 µl	90 µl

Freeze unused STET++ at -20°C until further use

## **Commented Protocol:**

### **1. Pick colony and grow in 3 ml LB + antibiotic over night at 37°C**

3 ml LB are fine if you prep 1.5 ml, use 4 ml if you want to prep 3 ml. Then you have still enough to inoculate a maxiprep, 500 µl for a glycerol stock or similar things...

### **2. Prepare the STET++ and switch on the heat block to 95 °C**

The STET++ should always be prepared fresh - and I always do. But the frozen rest works fine as well, so probably you can also prepare STET++ and freeze in aliquots.

### **3. Spin down 1.5 ml for 2 min at 8000 rpm (6000 G) store the rest at 4 °C**

Most other protocols recommend longer and higher centrifugation steps, but this step gives you a pellet that is easy to resuspend and would not fall off during the following steps. Do not pipet off the supernatant - open the eppi, discard all liquid inside and beat the eppi hard and several times upside-down on a piece of paper towel! The pellet will stand it and the liquid is efficient and fast removed.

### **4. Resuspend pellet in 300 µl STET++ buffer**

That is the most annoying step in the protocol. You have to resuspend each pellet with a fresh tip. Do not leave any pieces of the pellet undissolved, or the lyses will be incomplete.

### **5. Incubate at room temperature for 10 min**

That is the minimum time. I incubated for longer times up to one hour without any problems (for example because of many samples).

### **6. Incubate at 95 °C (in heat block) for 2 min (not longer!)**

That is really important. All protocols warn to not extend this step - so I never did.

### **7. Spin down for 15 min at 13000 rpm (18000 G)**

Some do spin shorter, but than you will not get a nice pellet easy to remove.

### **8. Discard pellet(!) with a toothpick a toothpick**

Works perfect, because the pellet sticks to it. But the first 1000 minipreps I picked the pellet with a yellow tip, which works with some training.

### **9. Add 500 µl Isoprop/NH4Ac and mix**

Just open all eppis, use one tip and pipet into all eppis, close them and shake the whole eppi-rack strongly.

### **10. Incubate for 20 min at room temperature**

This is also not time critical - The longer the better. But do not freeze it, otherwise NH4Ac may precipitate.

### **11. Spin down for 15 min at 13000 rpm (18000 G)**

Probably 10 min are sufficient. **You will not see a pellet!!!** Do not panic, everything is fine, the pellet is fine and transparent, it sticks to the wall and I never lost one - but it's invisible.

### **12. Discard supernatant and wash with 500 µl 70% Ethanol**

Like before, do not pipet off the supernatant - open the eppi, discard all liquid inside by turning it over and beat the eppi hard and several times upside down on a piece of paper towel! The pellet will stand it and the liquid is efficient and fast removed.

### **13. Spin down for 10 min at 13000 rpm (18000 G)**

Probably 5 min are sufficient.

### **14. Discard supernatant**

Like before, do not pipet off the supernatant - open the eppi, discard all liquid inside by turning it over and beat the eppi hard and several times upside down on a piece of paper towel! The pellet will stand it and the liquid is efficient and fast removed.

### **15. Spin down for 2 min at 13000 rpm (18000 G)**

Enough to spin down the remaining rest.

### **16. Pipett remaining supernatant of**

That is the most efficient way to remove the rest. You can also dry or speedvac, but it will take a long time and over dried DNA does not dissolve well.

### **17. Let pellet dry for 1 min**

If you are finished with pipetting, the DNA is dry enough. If you extend this to 15 min it would not be a problem as well.

### **18. Resuspend pellet in 30 µl demi water**

The DNA solution will be around 0.2 µg/µl. Allow some time to dissolve and pipet up and down a few times when you take out the DNA for the digest.

### **19. Digest 5 µl of DNA solution in a total volume of 30 µl (MasterMix !)**

Remember, this is a digest for testing the DNA, not for further cloning! So you need low amounts of restriction enzyme. I always use 0.3 µl enzyme in this assay, which is more than enough (if you use 0.1 µl all DNA will be digested after 5-10 min). Calculate how much samples you have, add 3 and make a MasterMix for that amount! I always incubate one hour - probably half an hour will do the job, too. Do not use expensive Enzymes. Sty I is a good choice. It will cut several times and the pattern of the bands will tell you if the clone is correct. **For more information on restriction digests, check the special protocols.**

### **20. Load 10 µl of the digest on agarose gel**

Remember: 0.7% for bands from 12.000 to 5.000 bp, 0.8% for 8.000 to 2.000 bp, 1.0% for 6.000 to 1.200 bp and 1.5% for 2.500 to 200 bp.

## **References and Comments:**

I got this protocol from Heinz Schaller in the ZMBH, Heidelberg. I did it many thousand times and never had any problems. Many people say that the DNA can not be digested by some restriction enzymes, but I never had any problems. I even used it for sequencing without problems, but then it's probably better to do a phenol/chloroform extraction first.

## **How to cite this page in publications:**

This document can be cited like this:

Untergasser, Andreas. "DNA-Miniprep. - Rapid boiling" *Untergasser's Lab*. Summer 2006.  
(include here the date when you accessed these page).

<[http://www.untergasser.de/lab/protocols/miniprep\\_rapid\\_boiling\\_v1\\_0.htm](http://www.untergasser.de/lab/protocols/miniprep_rapid_boiling_v1_0.htm)>.

## **Please Do Not Reprint This Article:**

This article is copyrighted. Please do not reproduce this article in whole or part, in any form, without [obtaining my written permission](#).

© by A. Untergasser

---

[Contact](#)

---

[Impressum & Disclaimer](#)